

How to Use the EKO Stethoscope

1. Click the purple 'Webside Connect' icon to submit a page request.



2. After full video activation by provider, user will be directed to use the EKO stethoscope.
3. User will press the middle button on the edge of the tablet to return to the main screen of the tablet.
4. Activate the EKO icon located on the main screen.
5. The provider will walk the bedside team through the process for stethoscope placement.
6. Post video encounter, bedside team will end the application, and close out of the session.

Best practice: User will connect the stethoscope to charger to ensure optimal charge at all times.